

## Product Datasheet

### Hoxd9 Antibody (orb1266604)

**Description**

Hoxd9 Antibody

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse

**Conjugation**

Unconjugated

**Tested**

WB

**Applications**
**Immunogen**

This Mouse Hoxd9 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 125-152 amino acids from the Central region of mouse Hoxd9.

**Target**

Hoxd9

**Preservatives**

Supplied in PBS with 0.09% (W/V) sodium azide.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

For WB starting dilution is: 1:1000

**Isotype**

Rabbit Ig

**Clonality**

Polyclonal

**MW**


35 kDa

**Uniprot ID**
[P28357](#)
**NCBI**
[P28357](#)
**Dilution Range**


For WB starting dilution is: 1:1000

**Expiration Date**

12 months from date of receipt.



Western Blot at 1:1000 dilution + mouse ...



Western blot analysis in 293 cell line I...