


Product Datasheet

SNTB1 Antibody (orb1266698)

Description	SNTB1 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	WB	Western blot analysis in NCI-H460 cell l...
Immunogen	This SNTB1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 243-271 amino acids from the Central region of human SNTB1.	
Target	SNTB1	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	Western blot analysis in mouse spleen ti...
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	58 kDa	
Uniprot ID	Q13884	
NCBI	Q13884	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	