
Product Datasheet

PITPNM1 Antibody (orb1266992)

Description

PITPNM1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

WB

Immunogen

This PITPNM1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 143-171 amino acids from the N-terminal region of human PITPNM1.

Target

PITPNM1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW



135 kDa

Uniprot ID
[O00562](#)
NCBI
[O00562](#)
Dilution Range


For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.

Western blot analysis in ZR-75-1 cell li...



Western blot analysis in mouse NIH-3T3 c...