

Product Datasheet

TOP3B Antibody (orb1267335)

Description	TOP3B Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This TOP3B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 826-854 amino acids from the C-terminal region of human TOP3B.
Target	TOP3B
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	97 kDa
Uniprot ID	O95985
NCBI	O95985
Dilution Range	For WB starting dilution is: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis in mouse bladder t...