

---

## Product Datasheet

### Map2k4 Antibody (orb1267944)

**Description**

Map2k4 Antibody

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Zebrafish

**Conjugation**

Unconjugated

**Tested Applications**

IHC-P, WB

**Immunogen**

This Mouse Map2k4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 259-287 amino acids from the C-terminal region of mouse Map2k4.

**Target**

Map2k4

**Preservatives**

Supplied in PBS with 0.09% (W/V) sodium azide.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

**Isotype**

Rabbit Ig

**Clonality**

Polyclonal

**MW**


44 kDa

**Uniprot ID**
[P47809](#)
**NCBI**
[P47809](#)
**Dilution Range**


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

**Expiration Date**


12 months from date of receipt.



Western blot analysis in mouse NIH-3T3 c...



Western blot analysis in Y79 cell line l...



Western blot analysis in zebra fish brai...