

Product Datasheet

RBM22 Antibody (orb1267958)

Description
RBM22 Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This RBM22 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 165-193 amino acids from the Central region of human RBM22.

Target

RBM22

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


47 kDa

Uniprot ID
Q9NW64
NCBI
Q9NW64
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of lysates from ze...



Western blot analysis in Jurkat cell lin...



Western blot analysis of RBM22 using rab...