

Product Datasheet

MINPP1 Antibody (orb1268102)

Description
MINPP1 Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications
Immunogen

This MINPP1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 370-398 amino acids from the C-terminal region of human MINPP1.

Target

MINPP1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

55 kDa


Uniprot ID
[Q9UNW1](#)
NCBI
[Q9UNW1](#)
Dilution Range

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.

NCI-H460



Western blot analysis in NCI-H460 cell I...



Flow cytometric analysis of NCI-H460 cel...