
Product Datasheet

Clk2 Antibody (orb1268180)

Description

Clk2 Antibody

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

This Mouse Clk2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 66-93 amino acids from the N-terminal region of mouse Clk2.

Target

Clk2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

60 kDa

Uniprot ID
[O35491](#)
NCBI
[O35491](#)
Dilution Range


For WB starting dilution is: 1:1000

Expiration Date


12 months from date of receipt.



Western Blot at 1:1000 dilution + mouse ...



Western blot analysis of lysate from mou...



Western blot analysis in mouse lung tiss...