


Product Datasheet

GABRB2 Antibody (orb1268241)

Description	GABRB2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This GABRB2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 346-374 amino acids from the C-terminal region of human GAB ₂ .
Target	GABRB2
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:500
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	59 kDa
Uniprot ID	P47870
NCBI	P47870
Dilution Range	For WB starting dilution is: 1:500
Expiration Date	12 months from date of receipt.



Western Blot at 1:500 dilution + HepG2 w...



Western blot analysis in A549 cell line ...