

Product Datasheet

Pak4 Antibody (orb1268295)

Description

Pak4 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, WB

Applications
Immunogen

This Mouse Pak4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 187-216 amino acids from the Central region of mouse Pak4.

Target

Pak4

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


65 kDa

Uniprot ID
[Q8BTW9](#)
NCBI
[Q8BTW9](#)
Dilution Range


For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis in MCF-7 cell line...



Western blot analysis in Neuro-2a cell line...



Flow cytometric analysis of Neuro-2a cells...