


Product Datasheet

PPM1G Antibody (orb1268336)

Description	PPM1G Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This PPM1G antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 64-92 amino acids from the N-terminal region of human PPM1G.	
Target	PPM1G	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	Western blot analysis in HepG2 cell line...
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	59 kDa	
Uniprot ID	O15355	
NCBI	O15355	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	