


Product Datasheet

NUMBL Antibody (orb1268571)

Description	NUMBL Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	This NUMBL antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 572-600 amino acids from the C-terminal region of human NUMBL.
Target	NUMBL
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	65 kDa
Uniprot ID	Q9Y6R0
NCBI	Q9Y6R0
Dilution Range	For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western blot analysis in 293 cell line I...



NUMBL Antibody immunohistochemistry anal...