



Product Datasheet

Pkmyt1 Antibody (orb1268572)

Description	Pkmyt1 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	WB	Western blot analysis in HL-60 cell line...
Immunogen	This Mouse Pkmyt1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 114-143 amino acids from the N-terminal region of mouse Pkmyt1.	
Target	Pkmyt1	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	Western blot analysis of Pkmyt1 using ra...
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	54 kDa	
Uniprot ID	Q9ESG9	
NCBI	Q9ESG9	Western blot analysis in mouse heart tis...
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	