




Product Datasheet

Pak7 Antibody (orb1268574)

Description	Pak7 Antibody	
Species/Host	Rabbit	 Western blot analysis in K562 cell line ...
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	
Immunogen	This Mouse Pak7 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 128-157 amino acids from the N-terminal region of mouse Pak7.	 Western blot analysis in mouse stomach t...
Target	Pak7	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	 Mouse Pak7 Antibody immunohistochemistry...
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	81 kDa	
Uniprot ID	Q8C015	
NCBI	Q8C015	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	