


Product Datasheet

Stk26 Antibody (orb1268576)

Description	Stk26 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	WB	Western blot analysis in human placenta ...
Immunogen	This Mouse Mst4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 271-299 amino acids from the Central region of mouse Mst4.	
Target	Stk26	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	Western blot analysis in mouse stomach t...
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	47 kDa	
Uniprot ID	Q99JT2	
NCBI	Q99JT2	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	