



Product Datasheet

Cytochrome P450 8B1 antibody (orb126858)

Description	Rabbit polyclonal antibody to Cytochrome P450 8B1 (also known as CYP12, CYPVIII B1) which is a liver ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9UNU6
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from 2...



Immunohistochemical analysis of formalin...