



---

## Product Datasheet

### TRAF1 Antibody (orb1268666)

<b>Description</b>	TRAF1 Antibody	
<b>Species/Host</b>	Rabbit	
<b>Reactivity</b>	Human, Mouse	
<b>Conjugation</b>	Unconjugated	
<b>Tested Applications</b>	WB	
<b>Immunogen</b>	This TRAF1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 55-83 amino acids from the N-terminal region of human TRAF1.	Western blot analysis in NCI-H460 cell l...
<b>Target</b>	TRAF1	
<b>Preservatives</b>	Supplied in PBS with 0.09% (W/V) sodium azide.	
<b>Form/Appearance</b>	Liquid	
<b>Concentration</b>	batch dependent	
<b>Storage</b>	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
<b>Note</b>	For research use only	Western blot analysis in mouse liver tis...
<b>Application notes</b>	For WB starting dilution is: 1:1000	
<b>Isotype</b>	Rabbit Ig	
<b>Clonality</b>	Polyclonal	
<b>MW</b>	46 kDa	
<b>Uniprot ID</b>	<a href="#">Q13077</a>	
<b>NCBI</b>	<a href="#">Q13077</a>	
<b>Dilution Range</b>	For WB starting dilution is: 1:1000	
<b>Expiration Date</b>	12 months from date of receipt.	