


Product Datasheet

Nek3 Antibody (orb1268681)

Description	Nek3 Antibody	
Species/Host	Rabbit	Western blot analysis in mouse liver tis...
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	
Immunogen	This Mouse Nek3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 446-474 amino acids from the C-terminal region of mouse Nek3.	
Target	Nek3	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	Mouse Nek3 Antibody immunohistochemistry...
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	57 kDa	
Uniprot ID	Q9R0A5	
NCBI	Q9R0A5	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	