

Product Datasheet

FCGR2A Antibody (orb1268811)

Description	FCGR2A Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This FCGR2A antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 256-285 amino acids from the C-terminal region of human FCGR2A.	Western blot analysis in NCI-H292 cell l...
Target	FCGR2A	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	35 kDa	
Uniprot ID	P12318	
NCBI	P12318	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	