


Product Datasheet

TLR7 Antibody (orb1268820)

Description	TLR7 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This TLR7 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 899-926 amino acids from the C-terminal region of human TLR7.
Target	TLR7
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	121 kDa
Uniprot ID	Q9NYK1
NCBI	Q9NYK1
Dilution Range	For WB starting dilution is: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis in HepG2 cell line...