
Product Datasheet

PITPNM2 Antibody (orb1268851)

Description

PITPNM2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

This PITPNM2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 639-668 amino acids from the Central region of human PITPNM2.

Target

PITPNM2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


149 kDa

Uniprot ID
[Q9BZ72](#)
NCBI
[Q9BZ72](#)
Dilution Range

For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis in 293 cell line
I...