


## Product Datasheet

# EPB41 (phospho-Tyr660/418) antibody (orb126904)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Rabbit Polyclonal to EPB41 (phospho-Tyr660/418).   |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Mouse   |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | WB   |
| <b>Immunogen</b>           | Internal   |
| <b>Concentration</b>       | 1mg/ml   |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks.<br>For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <b>P11171</b>  |
| <b>Dilution Range</b>      | WB: 1:500-1:3000, ELISA: 1:1000  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |



Western blot analysis of extracts from H...