

Product Datasheet

THAP8 Antibody (orb1269108)

Description	THAP8 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This THAP8 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 72-100 amino acids from the Central region of human THAP8.	Western blot analysis in NCI-H292 cell l...
Target	THAP8	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	30 kDa	
Uniprot ID	Q8NA92	
NCBI	Q8NA92	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	