

Product Datasheet

LHFPL1 Antibody (orb1269116)

Description	LHFPL1 Antibody	 <p>h.placenta 95 36 28 17</p>
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This LHFPL1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 20-49 amino acids from the N-terminal region of human LHFPL1.	Western blot analysis in human placenta ...
Target	LHFPL1	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	24 kDa	
Uniprot ID	Q86WI0	
NCBI	Q86WI0	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	