
Product Datasheet

CBLC Antibody (orb1269504)

Description

CBLC Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P, WB

Immunogen

This CBLC antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 424-453 amino acids from the C-terminal region of human CBLC.

Target

CBLC

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:50~100
For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


52 kDa

Uniprot ID
[Q9ULV8](#)
NCBI
[Q9ULV8](#)
Dilution Range


For IHC-P starting dilution is: 1:50~100
For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Formalin-fixed and paraffin-embedded hum...



Western blot analysis in K562 cell line

...