


Product Datasheet

SPATA7 Antibody (orb1269613)

Description	SPATA7 Antibody	 Western blot analysis in A549 cell line ...
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This SPATA7 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 374-401 amino acids from the C-terminal region of human SPATA7.	
Target	SPATA7	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	68 kDa	
Uniprot ID	Q9P0W8	
NCBI	Q9P0W8	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	