


Product Datasheet

FUT4 Antibody (orb1269812)


Description	FUT4 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	This FUT4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 495-524 amino acids from the C-terminal region of human FUT4.
Target	FUT4
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	59 kDa
Uniprot ID	P22083
NCBI	P22083
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western blot analysis in NCI-H460 cell I...



FUT4 Antibody immunohistochemistry analy...



Flow cytometric analysis of NCI-H460 cel...