


Product Datasheet

HOXA1/B1/D1 antibody (orb126984)

Description	Rabbit Polyclonal to HOXA1/B1/D1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	C-terminal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9GZZ0, P49639, P14653
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...