




---

## Product Datasheet

### ACAP1 Antibody (orb1270349)

<b>Description</b>	ACAP1 Antibody	
<b>Species/Host</b>	Rabbit	
<b>Reactivity</b>	Human, Mouse	
<b>Conjugation</b>	Unconjugated	
<b>Tested Applications</b>	IHC-P, WB	Western blot analysis in ZR-75-1 cell li...
<b>Immunogen</b>	This ACAP1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 96-125 amino acids from the N-terminal region of human ACAP1.	
<b>Target</b>	ACAP1	
<b>Preservatives</b>	Supplied in PBS with 0.09% (W/V) sodium azide.	
<b>Form/Appearance</b>	Liquid	Western blot analysis in mouse spleen ti...
<b>Concentration</b>	batch dependent	
<b>Storage</b>	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
<b>Note</b>	For research use only	
<b>Application notes</b>	For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100	ACAP1 Antibody immunohistochemistry anal...
<b>Isotype</b>	Rabbit Ig	
<b>Clonality</b>	Polyclonal	
<b>MW</b>	82 kDa	
<b>Uniprot ID</b>	<a href="#">Q15027</a>	
<b>NCBI</b>	<a href="#">Q15027</a>	
<b>Dilution Range</b>	For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100	
<b>Expiration Date</b>	12 months from date of receipt.	