


Product Datasheet

MLK1/2 (phospho-Thr312/266) antibody (orb127035)

Description	Rabbit Polyclonal to MLK1/2 (phospho-Thr312/266).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Internal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P80192, Q02779
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...