


Product Datasheet

RPH3AL Antibody (orb1270623)

Description	RPH3AL Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	Western Blot at 1:1000 dilution + A549 w...
Immunogen	This RPH3AL antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 208-237 amino acids from the C-terminal region of human RPH3AL.	
Target	RPH3AL	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	RPH3AL Antibody immunohistochemistry ana...
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	34 kDa	
Uniprot ID	Q9UNE2	
NCBI	Q9UNE2	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	