

Product Datasheet

SEL1L Antibody (orb1270665)

Description

SEL1L Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This SEL1L antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 156-182 amino acids from the Central region of human SEL1L.

Target

SEL1L

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000
For IHC-P starting dilution is: 1:50~100
For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


89 kDa

Uniprot ID
[Q9UBV2](#)
NCBI
[Q9UBV2](#)
Dilution Range


For WB starting dilution is: 1:1000
For IHC-P starting dilution is: 1:50~100
For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis in MCF-7 cell line...



SEL1L antibody immunohistochemistry anal...



Flow cytometric analysis of MCF-7 cells

...