

Product Datasheet

HFM1 Antibody (orb1270715)

Description

HFM1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This HFM1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1293-1321 amino acids from the C-terminal region of human HFM1.

Target

HFM1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


163 kDa

Uniprot ID
A2PYH4
NCBI
A2PYH4
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.

MDA-MB231
250


Western blot analysis
in MDA-MB231 cell ...



HFM1 antibody
immunohistochemistry
analy...



Flow cytometric
analysis of MDA-
MB231 ce...