

Product Datasheet

BCO1 Antibody (orb1270803)

Description	BCO1 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This BCMO1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 263-292 amino acids from the Central region of human BCMO1.	Western blot analysis in WiDr cell line ...
Target	BCO1	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	63 kDa	
Uniprot ID	Q9HAY6	
NCBI	Q9HAY6	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	