



Product Datasheet

GAS2L1 Antibody (orb1270901)

Description	GAS2L1 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, WB	
Immunogen	This GAS2L1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 618-646 amino acids from the C-terminal region of human GAS2L1.	Western blot analysis in K562 cell line ...
Target	GAS2L1	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	Flow cytometric analysis of K562 cells (...)
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	73 kDa	
Uniprot ID	Q99501	
NCBI	Q99501	
Dilution Range	For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	