




Product Datasheet

CCL21 Antibody (orb1270951)

Description	CCL21 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	
Immunogen	This CCL21 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 23-52 amino acids from the Central region of human CCL21.	Western blot analysis in 293 cell line I...
Target	CCL21	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	CCL21 Antibody immunohistochemistry anal...
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	15 kDa	
Uniprot ID	O00585	
NCBI	O00585	
Dilution Range	For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	