

Product Datasheet

LGALS9B Antibody (orb1270973)

Description

LGALS9B Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

This LGALS9B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 160-188 amino acids from the Central region of human LGALS9B.

Target

LGALS9B

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

40 kDa

Uniprot ID

[Q3B8N2](#)

NCBI


[Q3B8N2](#)

Dilution Range


For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis in 293 cell line l...



Western blot analysis in mouse NIH-3T3 c...



Flow cytometric analysis of 293 cells (r...