

## Product Datasheet

### ALOX15 Antibody (orb1271015)

## Description

ALOX15 Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested Applications

IF, IHC-P, WB

### Immunogen

This ALOX15 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 610-643 amino acids from the C-terminal region of human ALOX15.

### Target

ALOX15

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

75 kDa

### Uniprot ID

[P16050](#)

### NCBI


[P16050](#)

### Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50

### Expiration Date


12 months from date of receipt.



Western Blot  
at 1:2000  
dilution Lane  
1: ...



Western blot  
analysis in  
Hela (lane 1),  
...



Formalin-  
fixed and  
paraffin-  
embedded  
hum...