

Product Datasheet

Klrb1c Antibody (FITC) (orb1271111)

Description	Klrb1c Antibody (FITC)
Species/Host	Mouse
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Target	Klrb1c
Preservatives	Phosphate-buffered aqueous solution, $\leq 0.09\%$ Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	P27814
Expiration Date	12 months from date of receipt.