


Product Datasheet

Ly6g Antibody (PE/Cy7) (orb1271138)

Description	Ly6g Antibody (PE/Cy7)
Species/Host	Rat
Reactivity	Mouse
Conjugation	PE/Cy7
Tested Applications	FC
Target	Ly6g
Preservatives	Phosphate-buffered aqueous solution, $\leq 0.09\%$ Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	P35461
Expiration Date	12 months from date of receipt.



C57Bl/6
bone marrow
cells were
stained w...