

Product Datasheet

Il17a Antibody (FITC) (orb1271264)

Description	Il17a Antibody (FITC)
Species/Host	Rat
Reactivity	Mouse, Rat
Conjugation	FITC
Tested Applications	FC
Target	Il17a
Preservatives	Phosphate-buffered aqueous solution, $\leq 0.09\%$ Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	Q62386 , Q61453
Expiration Date	12 months from date of receipt.