



Product Datasheet

PAK5/6 (phospho-Ser602/Ser560) antibody (orb127127)


Description	Rabbit Polyclonal to PAK5/6 (phospho-Ser602/Ser560).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9P286, Q9NQU5
Dilution Range	WB: 1:500-1:3000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of PAK5/6(Phospho...



Western blot analysis of extracts from L...



Western blot analysis of extracts from H...