


Product Datasheet

PIGY antibody (orb127140)

Description	Rabbit Polyclonal to PIGY.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	N-terminal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q3MUY2
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...