

Product Datasheet

Ctla4 Antibody (PE/Cy7) (orb1271408)

Description	Ctla4 Antibody (PE/Cy7)
Species/Host	Hamster
Reactivity	Mouse
Conjugation	PE/Cy7
Tested Applications	FC
Target	Ctla4
Preservatives	Phosphate-buffered aqueous solution, $\leq 0.09\%$ Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Armenian Hamster IgG
Clonality	Monoclonal
Uniprot ID	P09793
Expiration Date	12 months from date of receipt.