


Product Datasheet

CD19 Antibody (PE) (orb1271496)

Description	CD19 Antibody (PE)
Species/Host	Mouse
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Target	CD19
Preservatives	Phosphate-buffered aqueous solution, ≤0.09% Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P15391
Expiration Date	12 months from date of receipt.



Human peripheral blood lymphocytes were ...