

Product Datasheet

FCGR3A Antibody (FITC) (orb1271681)

Description	FCGR3A Antibody (FITC)
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Target	FCGR3A
Preservatives	Phosphate-buffered aqueous solution, $\leq 0.09\%$ Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P08637
Expiration Date	12 months from date of receipt.