

Product Datasheet

IL2RA Antibody (PE/Cy7) (orb1271746)

Description	IL2RA Antibody (PE/Cy7)
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy7
Tested Applications	FC
Target	IL2RA
Preservatives	Phosphate-buffered aqueous solution, $\leq 0.09\%$ Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P01589
Expiration Date	12 months from date of receipt.