


Product Datasheet

CD4 Antibody (PerCP/Cy5.5) (orb1271919)

Description	CD4 Antibody (PerCP/Cy5.5)
Species/Host	Mouse
Reactivity	Human
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Target	CD4
Preservatives	Phosphate-buffered aqueous solution, ≤0.09% Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	P01730
Expiration Date	12 months from date of receipt.



Human peripheral blood lymphocytes were ...