



Product Datasheet

RNF144A antibody (orb127199)

Description	Rabbit Polyclonal to RNF144A.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Internal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P50876
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from C...



Western blot analysis of extracts from C...