



Product Datasheet

RSK1/2/3/4 (phospho-Ser221/227/S218/232) antibody (orb127205)

Description	Rabbit Polyclonal to RSK1/2/3/4 (phospho-Ser221/227/S218/232).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Internal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q15349 , Q15418 , P51812 , Q9UK32
Dilution Range	WB: 1:500-1:3000, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from H...



Western blot analysis of extracts from H...